

IHW#	Transfectant Name	Original Designation	Transfected	Antibiotic resista	Vector	Growth Medium	Comments
IHW03001	AMALA_A*24	AMALA	A*24	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03002	BM15_A*24	BM15	A*24	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03003	BM21_A*24	BM21	A*24	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03004	BSM_A*24	BSM	A*24	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03005	BTB_A*24	BTB	A*24	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03006	HOM2_A*24	HOM2	A*24	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03007	KAS001_A*24	KAS001	A*24	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03008	KOSE_A*24	KOSE	A*24	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03009	LWAGS_A*24	LWAGS	A*24	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03010	LZL_A*24	LZL	A*24	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03011	PF4015_A*24	PF4015	A*24	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03012	STEINLIN_A*24	STEINLIN	A*24	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03013	VAVY_A*24	VAVY	A*24	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03014	BM15_A*2	BM15	A*2	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03016	HOM2_A*2	HOM2	A*2	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03017	KAS116_A*2	KAS116	A*2	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03018	LO541265_A*2	LO541265	A*2	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03019	LWAGS_A*2	LWAGS	A*2	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03020	PF4015_A*2	PF4015	A*2	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03023	WT47_A*2	WT47	A*2	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03024	WTS1_A*2	WTS1	A*2	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03025	YAR_A*2	YAR	A*2	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03026	AMALA_A*3303	AMALA	A*3303	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03027	BM15_A*3303	BM15	A*3303	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03028	BM21_A*3303	BM21	A*3303	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03029	BOLETH_A*3303	BOLETH	A*3303	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03030	BSM_A*3303	BSM	A*3303	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03031	BTB_A*3303	BTB	A*3303	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03032	HOM2_A*3303	HOM2	A*3303	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03033	KAS001_A*3303	KAS001	A*3303	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03034	LO541265_A*3303	LO541265	A*3303	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03035	LWAGS_A*3303	LWAGS	A*3303	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03036	LZL_A*3303	LZL	A*3303	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03037	PF4015_A*3303	PF4015	A*3303	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03038	AMALA_B*4403	AMALA	B*4403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03039	BM15_B*4403	BM15	B*4403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03040	BM21_B*4403	BM21	B*4403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03042	BSM_B*4403	BSM	B*4403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03043	BTB_B*4403	BTB	B*4403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03044	HOM2_B*4403	HOM2	B*4403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03045	KAS001_B*4403	KAS001	B*4403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03046	LO541265_B*4403	LO541265	B*4403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03047	LWAGS_B*4403	LWAGS	B*4403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03048	LZL_B*4403	LZL	B*4403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03049	PF4015_B*4403	PF4015	B*4403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03050	AMALA_B*5201	AMALA	B*5201	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03051	BM15_B*5201	BM15	B*5201	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03052	BM21_B*5201	BM21	B*5201	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03053	BOLETH_B*5201	BOLETH	B*5201	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03054	BSM_B*5201	BSM	B*5201	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03055	BTB_B*5201	BTB	B*5201	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03056	HOM2_B*5201	HOM2	B*5201	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03057	KAS001_B*5201	KAS001	B*5201	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03058	LO541265_B*5201	LO541265	B*5201	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03059	LWAGS_B*5201	LWAGS	B*5201	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03060	LZL_B*5201	LZL	B*5201	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03061	PF4015_B*5201	PF4015	B*5201	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03062	AMALA_C*1403	AMALA	C*1403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03063	BM15_C*1403	BM15	C*1403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03064	BM21_C*1403	BM21	C*1403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03065	BOLETH_C*1403	BOLETH	C*1403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03066	BSM_C*1403	BSM	C*1403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03067	BTB_C*1403	BTB	C*1403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03068	HOM2_C*1403	HOM2	C*1403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03069	KAS001_C*1403	KAS001	C*1403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03070	LO541265_C*1403	LO541265	C*1403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03071	LWAGS_C*1403	LWAGS	C*1403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03072	LZL_C*1403	LZL	C*1403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03073	PF4015_C*1403	PF4015	C*1403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03074	AMALA_C*12	AMALA	C*12	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03075	BM15_C*12	BM15	C*12	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03076	BM21_C*12	BM21	C*12	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03078	BSM_C*12	BSM	C*12	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03079	BTB_C*12	BTB	C*12	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03080	HOM2_C*12	HOM2	C*12	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03081	KAS001_C*12	KAS001	C*12	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03082	LO541265_C*12	LO541265	C*12	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03083	LWAHS_C*12	LWAHS	C*12	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03084	LZL_C*12	LZL	C*12	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03085	PF401265_C*12	PF401265	C*12	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03086	.221_A*1101	.221	A*1101	Puromycin			
IHW03087	.221_A*24	.221	A*24	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03088	.221_A*0201	.221	A*0201	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03089	.221_A*2603	.221	A*2603	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03090	.221_A*3303	.221	A*3303	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03091	.221_B*4403	.221	B*4403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03092	.221_B*5201	.221	B*5201	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03093	.221_C*1202	.221	C*1202	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03094	.221_C*1403	.221	C*1403	Puromycin		use 0.6 ug/mL puromycin every 1-2 month for 3 days	
IHW03096	.221_C*0401	.221	C*0401	Puromycin		stable transfectant	
IHW03098	1331-8233_A* 2402	1331-8233	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03099	1331-8240_A* 2402	1331-8240	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03100	1332-8260_A* 2402	1332-8260	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03101	1332-8262_A* 2402	1332-8262	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03102	1333-8281_A* 2402	1333-8281	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03103	1333-8282_A* 2402	1333-8282	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03104	1341-8314_A* 2402	1341-8314	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03105	1341-8317_A* 2402	1341-8317	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03106	1344-8353_A* 2402	1344-8353	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03107	1344-8356_A* 2402	1344-8356	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03108	1346-8365_A* 2402	1346-8365	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03109	1346-8438_A* 2402	1346-8438	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03110	1347-8440_A* 2402	1347-8440	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03111	1347-8445_A* 2402	1347-8445	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03112	1349-8398_A* 2402	1349-8398	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03113	1349-8399_A* 2402	1349-8399	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03114	1362-8563_A* 2402	1362-8563	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03115	1362-8574_A* 2402	1362-8574	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03116	1408-1012_A* 2402	1408-1012	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03117	1408-1014_A* 2402	1408-1014	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03118	1413-1084_A* 2402	1413-1084	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03119	1413-1089_A* 2402	1413-1089	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03120	1416-1188_A* 2402	1416-1188	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03121	1416-1195_A* 2402	1416-1195	A* 2402	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03122	1331-8233_B*4403	1331-8233	B*4403	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03123	1331-8240_B*4403	1331-8240	B*4403	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03124	1332-8260_B*4403	1332-8260	B*4403	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03125	1332-8262_B*4403	1332-8262	B*4403	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03126	1333-8281_B*4403	1333-8281	B*4403	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03127	1333-8282_B*4403	1333-8282	B*4403	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03128	1341-8314_B*4403	1341-8314	B*4403	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	
IHW03129	1341-8317_B*4403	1341-8317	B*4403	Puromycin		use 0.7 ug/mL puromycin every 2 month for 5 days	

